

Möller Medical GmbH
Wasserkuppenstraße 29-31
36043 Fulda, Germany
Tel.: +49 (0) 661-9 41 95-0
Fax: +49 (0) 661-9 41 95-850
info@moeller-medical.com
www.moeller-medical.com

en

Instruction leaflet

30006443 / 2020-10 C



Instructions for Use TRAUMACEM™ V+ Injection Cannula

Manufacturer:
Möller Medical GmbH

Distributed by:



Device in scope

03.702.121S | TRAUMACEM™ V+ Injection
Cannula, for TFNA System

Description

The TRAUMACEM™ V+ Injection Cannula, for TFNA System is a sterile packed injection cannula designed for the augmentation of suitable DePuy Synthes implants.

Intended use

The TRAUMACEM™ V+ Injection Cannula, for TFNA System is intended for the augmentation of suitable TFN-ADVANCED® Proximal Femoral Nailing System Head Elements (TFNA Helical Blades and TFNA Screws) with PMMA¹-based bone cement in the proximal femur region. The injection cannula may only be used with the TRAUMACEM™ V+ Injectable Bone Cement and the TRAUMACEM™ V+ Syringe Kit.

¹ PMMA = Polymethyl Methacrylate

en

Refer to the corresponding instruction for use regarding indications, contraindications, compatibility, use, precautions, warnings and side effects of the bone cement used in conjunction with the TRAUMACEM™ V+ Injection Cannula. For information on compatibility with other devices or systems, consultation with a DePuy Synthes representative is recommended.

Indications

The indications for the augmentation procedure using the TRAUMACEM™ V+ Injection Cannula are defined by the PMMA-based bone cement and can be found in the respective instructions for use.

Contraindications

This device is not to be used in patients with coagulation disorders or infections. Further contraindications are defined by the PMMA-based bone cement and can be found in the respective instructions for use.

Patient Target Group

The patient target group is defined by the PMMA-based bone cement and can be found in the respective instructions for use.

Warnings and precautions

Since the device is used in technically complex surgery, the physician should be familiar with handling and use of the device, the other instrumentation and the procedure. It is also necessary to check whether the patient has allergies to materials (e.g. nickel) contained in the device. Should this be the case, the use of the device is at the

en

discretion of the attending physician.

The safe use of the device cannot be guaranteed if used in conjunction with MRI, which is associated with certain risks, including heating or migration of the device and artefacts on the MRI image.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established (only for European Union). Other countries should follow local regulations or procedures.

Side effects

The side effects are defined by the PMMA-based bone cement and can be found in the respective instructions for use.

Materials

The device contains the following materials:

- Stainless steel: Cannula including sleeve and plunger
- Polybutylene terephthalate: Handles and other plastic components

Device content

Each device contains a side-opening cannula (with Luer-lock) including a sleeve and a plunger.

Mode of use

For further information on the use of the TRAUMACEM™ V+ Injection Cannula, for TFNA System please refer to the DePuy Synthes surgical technique for TFN-ADVANCED® Proximal Femoral Nailing System.

en

Sterility

The device is supplied sterile. Before opening, check that the packaging is intact. Remove items from the packaging while maintaining aseptic conditions. Do not use if the packaging seal is broken or if the sterile barrier has been breached.

Single-use

The device is intended for single-use only and may not under any circumstances be reused. The reuse of a single-use device poses an infection risk to patients and users. Contamination of the device can lead to injury, illness or death of the patient. Cleaning, disinfection and sterilization can adversely affect important material and design properties and product function.

Storage

The device should be stored in a clean, dry environment and protected from direct sunlight. Do not use after the expiration date on the label.

Disposal

Disposal takes place in accordance with the hospital ordinance or the applicable laws. Due to the risk of infection and injury, care must be taken to avoid all contact with the tips or sharp edges of the instruments. This is especially important for contaminated instruments.

en

Symbols



Consult instructions for use



Batch code



Catalogue number



Keep dry



Keep away from sunlight



Manufacturer



Use by YYYY-MM-DD



Sterilized using ethylene oxide



Do not reuse



Do not resterilize



Do not use if the packaging is damaged



Only for USA: Caution: Federal law restricts this device to sale by or on the order of a physician.

For further information about symbols used please refer to our homepage: www.moeller-medical.com/glossary-symbols

en