


# Automatic Biopsy Gun & Coaxial Needles

Automatic disposable biopsy device for obtaining high quality histological samples

- Centimetre markings for precise determination of penetration depth
- 10 mm side notch
- Echogenic tip for clear visibility using ultrasound
- Clear colour-coding of the gauge sizes
- Can be used with one hand (ideal for prostate biopsies)
- 2 triggers (on the side and at the back)
- Choice of semi- or fully-automatic mode
- Strong spring
- Special low-noise damping for improved patient comfort
- Corresponding coaxial needles are to be found in the table below



Automatic Biopsy Gun ABG-2030

Biopsy needle size	Biopsy needle catalogue ref.	PU	Coaxial needle size	Coaxial needle catalogue ref.	PU
20G x 10 cm	ABG-2030-20-0100	5	19G x 5.1 cm	COX-2000-20-0051	10
20G x 15 cm	ABG-2030-20-0150	5	19G x 10.1 cm	COX-2000-20-0101	10
20G x 20 cm	ABG-2030-20-0200	5	19G x 15.1 cm	COX-2000-20-0151	10
18G x 10 cm	ABG-2030-18-0100	5	17G x 5.1 cm	COX-2000-18-0051	10
18G x 15 cm	ABG-2030-18-0150	5	17G x 10.1 cm	COX-2000-18-0101	10
18G x 20 cm	ABG-2030-18-0200	5	17G x 15.1 cm	COX-2000-18-0151	10
18G x 25 cm	ABG-2030-18-0250	5	17G x 20.1 cm	COX-2000-18-0201	10
16G x 10 cm	ABG-2030-16-0100	5	15G x 5.1 cm	COX-2000-16-0051	10
16G x 15 cm	ABG-2030-16-0150	5	15G x 10.1 cm	COX-2000-16-0101	10
16G x 20 cm	ABG-2030-16-0200	5	15G x 15.1 cm	COX-2000-16-0151	10
14G x 10 cm	ABG-2030-14-0100	5	13G x 5.1 cm	COX-2000-14-0051	10
14G x 15 cm	ABG-2030-14-0150	5	13G x 10.1 cm	COX-2000-14-0101	10
14G x 20 cm	ABG-2030-14-0200	5	13G x 15.1 cm	COX-2000-14-0151	10

Coaxial needle to guide the biopsy needles and to allow multiple sampling with a single puncture

COX-2000 is compatible with Automatic Biopsy Gun ABG-2030


- Centimetre markings and infinitely variable spring clips allow precise determination of the penetration depth
- Echogenic tip for clear visibility using ultrasound
- Hubs colour-coded to match the corresponding ABG disposable biopsy guns



## Automatic Biopsy Gun & Coaxial Needles

Automatic disposable biopsy device for obtaining high quality histological samples

- Centimetre markings for precise determination of penetration depth
- 20 mm side notch
- Echogenic tip for clear visibility using ultrasound
- Clear colour-coding of the gauge sizes
- Can be used with one hand (ideal for prostate biopsies)
- 2 triggers (on the side and at the back)
- Strong spring
- Special low-noise damping for improved patient comfort
- Fully automatic mode
- Corresponding coaxial needles are to be found in the table below



Automatic Biopsy Gun ABG-2040

Biopsy needle size	Biopsy needle catalogue ref.	PU	Coaxial needle size	Coaxial needle catalogue ref.	PU
20G x 10 cm	ABG-2040-20-0100	5	19G x 5.1 cm	COX-2000-20-0051	10
20G x 15 cm	ABG-2040-20-0150	5	19G x 10.1 cm	COX-2000-20-0101	10
20G x 20 cm	ABG-2040-20-0200	5	19G x 15.1 cm	COX-2000-20-0151	10
18G x 10 cm	ABG-2040-18-0100	5	17G x 5.1 cm	COX-2000-18-0051	10
18G x 15 cm	ABG-2040-18-0150	5	17G x 10.1 cm	COX-2000-18-0101	10
18G x 20 cm	ABG-2040-18-0200	5	17G x 15.1 cm	COX-2000-18-0151	10
18G x 25 cm	ABG-2040-18-0250	5	17G x 20.1 cm	COX-2000-18-0201	10
16G x 10 cm	ABG-2040-16-0100	5	15G x 5.1 cm	COX-2000-16-0051	10
16G x 15 cm	ABG-2040-16-0150	5	15G x 10.1 cm	COX-2000-16-0101	10
16G x 20 cm	ABG-2040-16-0200	5	15G x 15.1 cm	COX-2000-16-0151	10
14G x 10 cm	ABG-2040-14-0100	5	13G x 5.1 cm	COX-2000-14-0051	10
14G x 15 cm	ABG-2040-14-0150	5	13G x 10.1 cm	COX-2000-14-0101	10
14G x 20 cm	ABG-2040-14-0200	5	13G x 15.1 cm	COX-2000-14-0151	10

Coaxial needle to guide the biopsy needles and to allow multiple sampling with a single puncture

COX-2000 is compatible with Automatic Biopsy Gun ABG-2040

- Centimetre markings and infinitely variable spring clips allow precise determination of the penetration depth
- Echogenic tip for clear visibility using ultrasound
- Hubs colour-coded to match the corresponding ABG disposable biopsy guns

